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### ELECTRONIC BABY SCALE

#### MODEL 6401

Thank you for purchasing a Momert product. We hope you will be satisfied with our product throughout its service life. Your electronic baby scale is a sensitive measuring instrument, therefore please protect it from overload and from big dynamic loads! Never leave your baby unattended on the tray!

Please study the entire Operating Manual carefully before you start using the product. Keep the manual in a safe place for future reference. Make sure other people using the product are familiar with these instructions. Your new electronic baby scale has been designed for measuring your baby's weight and displaying his/her weight difference.

#### IMPORTANT SAFETY INSTRUCTIONS:

- Make sure you follow the marked polarity of batteries when inserting and replacing it.
- If the appliance will be out of use for a longer time, remove the batteries.
- Do not leave the unit unattended when in use.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
- Use the product only for the intended purpose specified in the Operating Manual.
- Use the unit only on a stable horizontal surface.
- Before cleaning, assembly, or disassembly of the cover, switch the unit off and remove the batteries!
- Do not use abrasive or chemically aggressive substances for cleaning the unit.
- Do not use the unit if it does not operate properly, or if it has been dropped, damaged or exposed to any liquids. Take the unit to an authorized service center for testing and repair.
- The unit is suitable for domestic use only. It is not designed for commercial operation.
- Keep the unit away from heat sources, such as radiators, ovens etc. Do not expose the unit to direct sunlight or moisture.
- Do not touch the unit with wet or damp hands.

- The product is a fragile measuring device. Prevent it from crashing or falling.
- Do not repair the appliance yourself. Contact an authorized service center.

#### FAILURE TO FOLLOW THE MANUFACTURER'S INSTRUCTIONS MAY LEAD TO REFUSAL OF WARRANTY REPAIR!

#### TECHNICAL SPECIFICATION:

- Automatic hold function
- Precision measuring system
- Weight difference display function
- Tare function
- Auto power off
- Low battery indicator
- Overload indicator

Technical parameters	
<b>Voltage</b>	2 pcs 3V CR2032 lithium batteries
<b>Weighing range</b>	20 kg
<b>Weighing increments</b>	5 g
<b>Weighing units</b>	Kg, Lb, St
<b>Working/Storage temperature</b>	20-25°C/15-40°C
<b>Display size</b>	LCD diameter= 8,0 cm , 3,15 inch

#### PRODUCT DESCRIPTION:

1. Weighing surface
2. ON/OFF button
3. UNIT button
4. TARE button
5. Display

#### OPERATING BUTTONS:

- The scale can be switched on and off with the **ON/OFF button**. (Push the button once.)
- With the **TARE button** you can use the tare function. "0.00" appears on the display 2 seconds after pressing the button.
- **UNIT button** is suitable to choose the required unit of measurement, the actual unit of

measurement can be seen on the right side of the display.

#### INSTRUCTION MANUAL:

1. The scale includes replaceable lithium batteries. Please remove the protective film of the battery before use (Fig. 1).

Fig. 1

2. Please set the unit of measurement (kg-lb-st) using the UNIT button.

Fig. 2

3. Place the scale on a hard, flat surface. (Do not use it on soft surface)

4. By pressing the ON/OFF button once the scale turns on. Please wait until "0,00" kg is displayed (Fig. 2).

5. Put the baby on the scale. The measured weight will be displayed in a few seconds, and it will be saved automatically in the memory. When there is an "M" on the display, then there is stored information in the memory. Don't leave the baby unattended on the scale.

6. By pressing the ON/OFF button once the scale turns off. Scale turns off automatically after 30 seconds if there is no measurement made.

After the next measurement the display shows in 2 seconds switching the measured weight and the difference to the previous measurement. Last measured weight is stored automatically in the memory.

#### MEMORY:

When there is an "M" on the display if the scale is turned on, it means the last measured weight is stored in the memory to be able to display the weight difference. The scale is suitable for measuring the quantity of the consumed milk, if you weight the baby before and after breast-feeding. If you don't need the weight difference to be shown, then by pressing the Tare button for 3 seconds the memory can be deleted. "M" will disappear from the display. If you take out the batteries, then the memory will also

be deleted.

#### NOTICE:

Temperature can effect the measuring data. An accurate measuring can be made at room temperature (20 - 25 °C). If you not only use the scale for measuring, but you put a blanket or a pillow also in it, then you have to set zero with Tare button before measuring the baby. If the weight of these accessories are less then 5 kg, then the maximum measured weight can be 20 kg. If the weight of them is more then 5 kg, then the maximum measured weight can be less then 20 kg. You have to discount the added weight of the 20 kg and this will be the maximum measurable weight.

#### OTHER INFORMATION ON THE DISPLAY:

If the "Lo" warning is displayed after switching on (Fig. 3), or there is nothing on the display, the battery is dead and has to be replaced.

Fig. 3

If the sign „ERR" appears on the display (Fig. 4.), it means that the scale is overloaded, in this case the weight must be reduced ( an overload can cause sensor faults).

Fig. 4

If the message „OUT2" appears on the display (Fig. 5.), it means sensor fault ( after overloading or after a sudden, big dynamic load ). In this case, we have no warranty liability for the product.

Fig. 5

#### CHANGE BATTERY:

Please open the battery holder at the backside of the scale and carefully take out the battery by pressing the small silver colour holding plate to get the battery out of the battery holder.

Replace the low battery to a new battery. The low battery should be thrown to the appropriate place.

#### BATTERIES ARE NOT INCLUDED IN WARRANTY! CLEANING AND MAINTENANCE:

The surface of the unit may only be cleaned using a wet cloth. Do not use any detergents or hard objects as they may damage the unit's surface!

**Never clean the unit under running water, do not rinse it, and do not submerge it in water!**

#### REPAIR AND MAINTENANCE:

Any extensive maintenance or repair that requires access to the inner parts of the product must be performed by an expert service centre.

#### ENVIRONMENTAL PROTECTION:

- Packaging materials and obsolete appliances should be recycled.
- The transport box may be disposed of as sorted waste.
- Polyethylene bags shall be handed over for recycling.

#### Appliance recycling at the end of its service life:

A symbol on the product or its packaging indicates that this product should not go into household waste. It must be taken to the collection point of an electric and equipment recycling facility. By making sure this product is disposed of properly, you will help prevent the negative effects on the environment and human health that would otherwise result from inappropriate disposal of this product. You can learn more about recycling this product from your local authorities, a household waste disposal service or in the shop where you bought this product.

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